



ANALIZA BIAŁEK



100003038

## Laboratory of Analysis

**Client:** Marika Mej  
**Order number:** 100003038  
**Name:** NASSA LABS OXANDROLONE  
**Number of samples:** 1 pill

**Batch number:** LX12/22  
**Recive date:** 2023-01-13  
**Start date:** 2023-01-19  
**End date:** 2023-01-19

## Sample identification

Substance	Method of taking a sample	Form
Oxandrolone	Sample provided by customer	White pills in a blister pack

## Results

Substance	Sample name	Content	Purity
Oxandrolone	NASSA LABS OXANDROLONE 10MG	9,86mg/tab	-----

## Research comment

The test sample was applied in the amount of 6µl to the BioBasic 8 Thermo Scientific column. Weight of the whole pill is 204,39mg and weighted portion is 173,87mg. During the course of the analysis with the method a linear gradient in the solutions of H2O and Acetonitrile with the addition of 0.1% trifluoroacetic acid, a peak of the sought substance appeared. Its content, when compared to the reference chromatogram, is 9,86mg/tab. Height the measurement error is + - 0,05mg. The actual state of the sample sent before the test is presented on the photo below

## Photo of the test sample



## Password

617DXCXA

**Analiza Białek sp. z o.o.**  
**KRS:** 0000966168  
**NIP:** 8992922970  
**REGON:** 521751875

**Headquarters and Laboratory:**  
ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**  
Analysis and research specialist  
mgr Arkadiusz Zajac




ANALIZA BIAŁEK



100003038

## Chromatogram



## Method

Time	H2O + 0,1% TFA	ACN + 0,1% TFA
0min	60%	40%
2min	60%	40%
27min	5%	95%
30min	5%	95%

## Password

617DXCXA

**Analiza Białek sp. z o.o.**

**KRS:** 0000966168

**NIP:** 8992922970

**REGON:** 521751875

**Headquarters and Laboratory:**

ul. Duńska 9, 54-427 Wrocław  
Budynek DELTA

**Lab technician:**

Analysis and research specialist  
mgr Arkadiusz Zajęc